

Algorithms integration

- 1 Definitions
- 2 Protocol
- 3 Test server and client
 - 3.1 Test server
 - 3.1.1 Running test server
 - 3.1.2 Behavior
 - 3.2 Test client
 - 3.2.1 Running client server
 - 3.2.2 Behavior
- 4 Requests
 - 4.1 Initialization
 - 4.1.1 Fields description
 - 4.2 Status update
 - 4.2.1 Fields description
- 5 Example request bodies
 - 5.1 Initialization
 - 5.2 Status update - error
 - 5.3 Status update - in_progress
 - 5.4 Status update - completed
 - 5.5 Request for tile

Definitions

CC.AI - CancerCenter.ai

TPA - Third-Party Algorithm

Protocol

Communication is done by HTTP requests. For one algorithm run there are following steps:

1. CC.AI starts by sending by sending Initialization request to TPA URL. TPA should respond with status code 200 and run algorithm on image using data from Initialization request.
2. The algorithm is in progress. TPA sends Status update requests that inform CC.AI about current progress of processing. The time between subsequent Status update requests (and Initialization and first Status update) should not exceed one minute - otherwise CC.AI will consider the algorithm run as crashed.
3. The algorithm returns the result. It should send the Status update request with status completed or error, indicating the result. After that, all Status update requests will be ignored.

Test server and client

Test server and client are provided to help with implementing protocol on TPA side - test-server imitates CC.AI and test-client is an example implementation of TPA. They are both written using Python. In order to run the server and the client user should install dependencies provided in requirements.txt .

Test server [↗](#)

Running test server [↗](#)

`test-server` is activated by running command `python main.py` inside `\server` directory. It supports following command-line arguments:

| Argument | Mandatory | Description | Default value |
|---------------------------|-----------|--|---|
| <code>--image_path</code> | Yes | Path to the image to process. It will be accessible on <code>tiles_url</code> from <code>Initialization</code> request | -- |
| <code>--tpa_url</code> | No | URL to which <code>test-server</code> will send <code>Initialization</code> request | <code>http://127.0.0.1:12346/run_algorithm</code> |
| <code>--auth</code> | No | If specified, its value will be set as <code>Authorization</code> header value while sending <code>Initialization</code> request | -- |
| <code>--port</code> | No | Port number which <code>test-server</code> will use | 12345 |
| <code>--host_url</code> | No | This URL will be used to create <code>return_url</code> field in <code>Initialization</code> request | <code>http://127.0.0.1</code> |

Example activation: `python main.py --image_path path/to/data/1.svs --port 54321`

Behavior [↗](#)

Upon activation it will send `Initialization` request to `--tpa_url` and start listening for `Status update` requests. For each of those requests it will print on command-line information that it received. When `Status update` with `status` field set to `completed` or `error` is received, `test-server` will print the result and finish. The server will also provide tiles from image specified by `image_path` argument by responding on endpoint specified in `tiles_url` field.

Test client [↗](#)

Running client server [↗](#)

`test-client` is activated by running command `python main.py` inside `\client` directory. It supports following command-line arguments:

| Argument | Mandatory | Description | Default value |
|---------------------|-----------|--|---------------|
| <code>--port</code> | No | Port number which <code>test-client</code> use | 12346 |

Example activation: `python main.py --port 54321`

Behavior [↗](#)

When ran, client will listen for `Initialization` requests, upon receiving one it will start processing it. It will download tiles from image, while sending `Status update` requests with progress updates in meantime. After that it will either randomly send `Status update` with status `error` or notify `CC.AI` about further progress and end processing with `Status update` request that has status `completed`, and sample result set.

Requests

Initialization

Initialization is a POST request with following body:

```
1 {
2   "image": {
3     "levels": int,
4     "width": int,
5     "height": int,
6     "tile_size": int,
7     "tiles_url": string,
8     "objective_magnification": float,
9     "microns_per_pixel": float
10  },
11  "return_url": string,
12  "id": string
13 }
```

Fields description

| Field | Description |
|-------------------------------|--|
| id | Unique ID assigned to single algorithm run |
| return_url | URL to which TPA should send Status update requests |
| image | Object that contains information about image on which the algorithm should be ran. The image is available to download as a tile pyramid. The file structure of the pyramid is described below, but for better understanding, one can search for information about Deep Zoom Image format (which we actually use under the hood). |
| image.levels | Number of zoom levels of image. Level <code>image.levels - 1</code> corresponds to original image, with size <code>image.height</code> x <code>image.width</code> . <code>image.levels - 2</code> corresponds to image downsampled 2 times, with size <code>image.height / 2</code> x <code>image.width / 2</code> , level <code>image.levels - 3</code> to image downsampled 4 times, with size <code>image.height / 4</code> x <code>image.width / 4</code> and so on. Level <code>0</code> corresponds to image downsampled enough times to have size 1x1 |
| image.width | Width of the image |
| image.height | Height of the image |
| image.tile_size | Width and height of the tiles returned by sending requests to <code>image.tiles_url</code> . Last tiles of either dimension can be smaller than this value |
| image.tiles_url | URL of form <code>https://example.com/path/{level}/{x}_{y}.jpeg</code> , that can be used to download image tiles. TPA should substitute desired values for <code>{level}</code> , <code>{x}</code> and <code>{y}</code> . <code>x</code> and <code>y</code> are tile coordinates meaning that call with values <code>level = image.levels - 2</code> , <code>x = 3</code> and <code>y = 4</code> will return tile starting at point(<code>x = 3 * image.tile_size</code> , <code>y = 4 * image.tile_size</code>) from image downsampled 2 times |
| image.objective_magnification | Objective magnification that was used to scan image (typically 10, 20 or 40) |
| image.microns_per_pixel | Microns per pixel of the scan |

Status update [↗](#)

Status update is a `POST` request with following body:

```
1 {
2   "status": str
3   "progress": int
4   "error": str
5   "result": object
6 }
```

Fields description [↗](#)

| Field | When mandatory | Description |
|----------|---|--|
| status | Always | Status update that <code>TPA</code> wants to inform about. Possible values are: <code>in_progress</code> , <code>error</code> and <code>completed</code> |
| progress | When status is <code>in_progress</code> | Progress of algorithm run. Number in range 0 - 100 |
| error | When status is <code>error</code> | Reason of failed processing |
| result | When status is <code>completed</code> | Object containing result of algorithm |

Example request bodies [↗](#)

Initialization [↗](#)

```
1 {
2   "image": {
3     "levels": 17,
4     "width": 41381,
5     "height": 52329,
6     "tile_size": 512,
7     "tiles_url": "http://127.0.0.1:54321/slide_files/{level}/{x}_{y}.jpeg",
8     "objective_magnification": 40,
9     "microns_per_pixel": 0.2491
10  },
11  "return_url": "http://127.0.0.1:54321/integrations/algorithm/5fb22a48-2e55-4e3f-8d91-1d5c32b67f38/status/",
12  "id": "5fb22a48-2e55-4e3f-8d91-1d5c32b67f38"
13 }
```

Status update - `error` [↗](#)

```
1 {
2   "status": "error",
3   "error": "Random error occured"
4 }
```

Status update - `in_progress` [↗](#)

```
1 {  
2   "status": "in_progress",  
3   "progress": 30  
4 }
```

Status update - completed [↗](#)

```
1 {  
2   "status": "completed",  
3   "result": {  
4     "regions_of_interest": [  
5       {  
6         "x": 20690,  
7         "y": 26164,  
8         "w": 100,  
9         "h": 100,  
10        "label": "Label 1",  
11        "score": 0.97  
12      }  
13    ],  
14    "other": "anything"  
15  }  
16 }
```

Request for tile [↗](#)

http://127.0.0.1:12345/slide_files/9/0_0.jpeg